Automated M² Measurement System





KEY FEATURES

1. LARGE APERTURES

The only M^2 system on the market equipped with a complete set of 50mm (2") optics. Also, the sensor is 11.3 x 11.3mm

2. SIMPLE ALIGNMENT

Two beam-steering mirrors are included for quick and easy alignment of your laser into the system. The internal mirrors are factory-aligned and the pre-set height also simplify the alignment

3. COMPACT

The low-profile ingenious mechanics make it easy to fit the device on any optical table

4. ISO COMPLIANT

The calculations are fully compliant to the ISO 11146 and 13694 standards

5. FAST ACQUISITION

Make a complete, ISO-compliant measurement in only 20 seconds with the ROI feature and in less than a minute with full-frame acquisition


6. FLEXIBLE & INTUITIVE SOFTWARE

In the easy-to-navigate software, both automatic and manual settings are available, so data points can be added or removed even after an automatic scan is completed

USER INTERFACE



Enter measurement parameters in the $\ensuremath{\mathsf{M}}^2$ Setup tab.



View and save results with the comprehensive $\ensuremath{\mathsf{M}}^2$ Results tab.

SEE ALSO

ACCESSORIES FOR BEAM DIAGNOSTICS LIST OF REGULAR ACCESSORIES

190 206

SPECIFICATIONS

BEAMAGE-M2

11.3 x 11.3 mm sensor

SENSOR TECHNOLOGY Beamage-4M included Ø 48 mm optics **EFFECTIVE APERTURE**

MEASUREMENT CAPABILITY

System Wavelength Range

3 Flip-mount attenuators for 8 levels of attenuation: no attenuation, ND0.5, ND1, ND2, ND1.5, ND2.5, ND3, ND3.5 Attenuation Range

Beam Diameter Range ^a 55 μm to 11.3/3 mm

Translation Stage

Mechanical Travel Range 200 mm Effective Optical Path Range 400 mm

Lens Focal Length 3 AR-coated lenses included: 200 mm, 250 mm and 300 mm

Typical M2 Accuracy b ±5% Typical M2 Repeatability b ±2%

Applicable Light Sources CW and pulsed

45 sec with full-frame acquisition Typical Measurement Time

DAMAGE THRESHOLDS °

Maximum Average Power 1 W with ND filter

CW: 10 W/cm²; Pulsed: 300 µJ/cm² Maximum Density (1064 nm)

PHYSICAL CHARACTERISTICS

Dimensions

Main Enclosure 357 mm (L) x 165 mm (W) x 135 mm (H) 602 mm (L) x 193 mm (W) x 172 mm (H) Total (including external mirrors)

Optical Axis Height 86 mm Weight

Power Supply 48V DC, 1.25A out

SOFTWARE

Displays 2D, 3D, XY, Beam Tracking and M²

Beam Diameter Definitions D4 σ

1/e² along crosshairs (13.5%) FWHM along crosshairs (50%)

Custom (%)

Beam Quality Definitions Laser beam quality M^2 : M_x^2 , M_y^2 (ISO compliant)

Beam Propagation Factor: BPP, BPP, Width at waist: W., W.

Waist location and offset: Z_{v} , Z_{v} , ΔZ Divergence angle: θ_{v} , θ_{v} Rayleigh length: Z_{Rx} , Z_{Ry}

Astigmatism

Printing and Reports Full report in print-ready format

ORDERING INFORMATION

Product Name Beamage-M2 Call Product Number

Specifications are subject to change without notice

Specifications in the table above are for the use with a Beamage-4M beam profiler (included in the Beamage-M2 kit)


- a. At the Beamage sensor
- Depending on the beam quality and optical configuration
- c. With ND4 filter at the Beamage

INTUITIVE SOFTWARE INTERFACE

Select which set of Rayleigh range boundaries to display on the graph: X, Y or both


M2 Results tab:

View and understand all the measured parameters quickly, for both the initial laser beam and the beam inside the BEAMAGE-M2 system



Switch to 3D or 2D displays to see each of the measured profiles

M2 Setup tab: Control your acquisition parameters



AUTOMATED MEASUREMENTS

Inside the BEAMAGE-M2, a computer-controlled motorized rail allows precise positioning of two folding mirrors, which in turn allow a 400 mm beam path difference. At each position of the translation stage, a beam profile is acquired and the beam diameter is measured. The automation of the translation stage allowed by the software is the key to a fast measurement.

AUTOMATED & EASY TO USE, YET POWERFUL & FLEXIBLE



ISO-COMPLIANCE MADE SIMPLE

With the "ISO SCAN" button, the software automatically defines new parameters for a more precise M^2 measurement. The "ISO SCAN" data set complies with the ISO-11146 M^2 measurement standard, being spread between $-3Z_R$ and $+3Z_R$.

The automatic settings are updated after each calculation, considering the values of Z_0 and Z_R from the latest measurement.

By default, the results graph always shows the calculated positions of the first three Rayleigh distances on each side of the waist. The "X" and "Y" button toggles them on or off.



FULL CONTROL ON YOUR DATA

During an M² scan, each of the measured profiles is saved and the flexible software gives you complete control on your acquired data.

- View each acquired profile in 2D display or 3D display.
- Add measurement points to a data set at the position of your choice with the "ADD" button.
- Remove unwanted profiles from your data set & recalculate the measurements.
- Change the beam diameter definition and the crosshair mode.



PRACTICAL ALIGNMENT TOOL

Each BEAMAGE-M2 system includes an alignment tube that helps you set up the system faster. Simply use the two alignment mirrors to center your laser beam onto both irises, and you will be ready to start measuring in no time!

The fluorescent material around the pinholes also helps to align beams that are in the NIR range without having to use an IR viewer.



FAST ATTENUATION

Add or remove attenuation with the flick of a finger. The software adjusts the exposure time at each frame during an acquisition, and advises the user on the required attenuation.